

# **Universiti Malaysia Perlis**

#### INVENTORS

DR. LIU WEI WEN PROF'DR. UDA HASHIM HAARINDRA PRASAD VEERADASAN PERUMAL

#### CONTACT DETAILS

institute of Nano Electronic Engineering (INEE) Universiti Malaysia Perlis 01000 Kangar, Perlis e-mail : wwliu@unimap.edu.my

# NANO DNA LAB-ON-CHIP MEDICAL DIAGNOSTIC SYSTEM FOR MULTI DETECTION OF LEPTOSPIROSIS



AVSIA PERLIS

## PRODUCT DESCRIPTION

The Nano DNA Lab-On-Chip is a medical diagnostic system with ultra-high sensitive and selective which is developed using carbon nanotubes/nanowires based, microfluidic-integrated biosensor for leptospirosis detection. Rapid detection at concentration as low as 1 pico molar and reusable for 5 times can be achieved by this device. Therefore, this device has huge potential in diagnostic not only for leptospirosis but other tropical diseases as well.

#### **INVENTION ADVANTAGES**

- · Biocompatible.
- · Small size (as big as one ten cent).
- Functionalized carbon nanotubes: binding to a broad variety of target molecule, good in affinity and specificity, easy and low cost to produce.
- Microfluidic: compact size, single cell handling and analysis, and micro droplets.

## **NOVELTIES**

- · Low cost fabrication of Nano DNA Lab-On-Chip.
- Label-free detection in static and real-time.
- Detection of ultra low concentration and high reusability.
- Rapid detection of leptospirosis.

## **COMMERCIAL POTENTIAL**

- · Highly potential for commercializing in clinical sample testing.
- Hand-held electronic device for rapid and selective detection.

#### APPLICATIONS

- A hand-held diagnostic device for medical practitioners for tropical diseases detection.
- A home life-care diagnostic appliance.






Figure 2: Functionalized carbon nanotubes



Figure 3: Comparison between DNA Lab-On Chip and ten cent sizes.



Figure 4: PDMS microfluidics.




Figure 5: pH measurement.




Figure 6: Measurement of glucose concentration.




Figure 7: Measurement of leptospirosis DNA concentration.





